

Shanghai Omni Laser Skinology CO.,LTD



EMT Muscle Sculpt Slimming System



ISO
13485



Health
Canada




APOSTILE

Aesthetic and Medical Equipment Manufacturer

EMT MUSCLE SCULPT SLIMMING SYSTEM

How It Works

RD-EMT muscle sculpt series uses High-intensity Focused Electromagnetic Teslasculpt technology, directly stimulates motor neurons, so that the body's muscles continue to expand and contract (this kind of contraction cannot be achieved by your usual exercise or fit exercise). The energy pulse of 30 minutes treatment can stimulate 30000 strong muscle contractions, which help the fat cells to metabolize and decompose vigorously. At the same time, with muscle strengthening, it brings new technological experience for body shaping.



Application

RD-EMT series is able to induce short bursts of powerful muscle contractions. The super-muscle exercises can achieve a dual effect. It not only deeply reshapes the internal structure of the muscle to increase the density and volume of the muscle, but also can quickly activate the movement of fat cells to achieve its supramaximal lipolysis. It is indicated to be used for:



01. ≡


Improvement of abdominal tone, strengthening of the abdominal muscles, development of firmer abdomen.

02. ≡

Strengthening, toning and firming of buttocks, thighs and calves.

03. ≡

Improvement of muscle tone and firmness, for strengthening muscles in arms.



Certificate



Advantage

Two Years' R&D Research

Different solutions for men and women muscle shaping. Safe and intelligent power supply utilized.

Five Exercise Modes

Simulate real fitness rhythm, effective exercise intensity and interval combination. It's like working out with pamel.

- ★ Hiit
- ★ Hypertrophy
- ★ Strength
- ★ Combo 1
- ★ Combo 2

Qualified Cooper Strip Coils

Omni Laser increases the number of turns of the coil, reduces the gap between the high energy. High energy use of cooper strip coils increases the number of turns of the coil, reduces the gap between the coil and the coil, and more effective increases the inductance, the measured inductance reaches 30uH, which can reach deeper muscle contractions. Long discharge time: while the inductance is increased, it can ensure that the output energy is more powerful, which makes the treatment more comfortable. Relatively low heat: the better matching air duct design allows long working time with high frequency and high energy, and continuous work can be achieved.

Omni Laser



Others






Intelligent Software

Simple mode and expert for your choice

Simple mode for beginners, frequency and intensity can be set by hand.

Expert mode for professional operator

Different work modes and different steps can be set by hand, and all operation steps can be saved in the programming, which can be uploaded automatically next time. 500 programming exercises can be saved in the system. Easy to manage.



FAQ

Q1: How does Muscle Stimulator build muscle and burn fat?

High Intensity Focused Electromagnetic enables muscle treatment to reach 20,000 "super muscle exercises" at one time, stimulates fat cell super-fast metabolism.

Strong fat decomposition, and promotes muscle strengthening and growth, muscle gain+ fat dissolving, in one step.


Q2: How many times is a course of treatment? How often do you do it at intervals?

Usually 6 times a course of treatment, 2 times a week, 30 minutes each time (to avoid menstrual period) after a course of treatment 19% reduction in fat, 16% increase in Muscle Mass.

Specification




Handles Choices




Model No.	RD-EMT600	RD-EMT700
Intended Use	1. Improvement of abdominal tone, strengthening of the abdominal muscles, development of firmer abdomen; 2. Strengthening, toning and firming of buttocks, thighs and calves; 3. Improvement of muscle tone and firmness, for strengthening muscles in arms.	
Principle of Action	Initiating action potential of nerves results in muscles contraction	
Display	15.6" colorful touch screen	12.1" colorful touch screen
Type of Energy	Magnetic field	
Number of Outputs	2 or 4	
Number of Magnetic Coils in Applicator	1	
Frequency	1-50HZ	
Therapy Time	30mins	


Before&After






Shanghai Omni Laser Skinology CO.,LTD

 Building 12-A-6F, Liando U Valley Songjiang International Enterprise Port Park, No. 1698 Minyi Road, Songjiang District, Shanghai, China

 Tel: +86-21-54847192 54847783

 info@omni-laser.com

 www.omni-lasers.com

 [www.youtube.com/
@OmniLaserSkinology](http://www.youtube.com/@OmniLaserSkinology)

