

# Shanghai Omni Laser Skinology CO.,LTD



## CET Body Health System





ISO 13485 **MDSAP** MEDICAL DEVICE SINGLE AUDIT PROGRAM **FDA**  Health Canada  **invimö** **CE** **PS E** **APOSTILE**

Aesthetic and Medical Equipment Manufacturer

# CET BODY HEALTH SYSTEM

## How It Works

The RD-CET device stimulates healthy cell function and revitalization through radio frequency therapy. It works at a specific frequency of 448kHz, CET body health system based on the medical studies described the Proionic System that stimulates intra and extra cellular ion exchange restoring physiological cell electrical activity, which enables the correct functionality of the cells, which causes thermal and sub-thermal effects, achieving noticeable results from the first session.




## Certificate



Face&body rejuvenation



Skin tightening&reduced orange peel tissue



Body shaping&slimming



Joint rehabilitation




Anti-hair loss



Improve dysmenorrhea&uterine cold



Fade gestation wrinkles and scars



Boost immunity



Reduce visceral fat




Eliminate inflammation and prevent disease



# Application



## Capacitive Handpiece(CAP)

- Skin and superficial tissue
- Improves orange peel tissue and relieves muscle pain
- Body shaping
- Lift firming
- Postoperative care, shorten the treatment time early wound recovery, shorten the period of internal bleeding edema



## Resistive Handpiece(RES)

- Deep heating
- Promote blood circulation throughout the body
- Promote visceral fat burning and so on
- Ameliorate cold symptoms
- Relieve inflammatory pain



## Advantages

### Pain-Free and Comfortable

Unlike some other body slimming treatments, using the 448khz CET machine is generally pain-free and comfortable. Many people find the warm sensation of the treatment to be quite relaxing.

# Non-invasive Body Contouring

The 448khz CET machine offers non-invasive body contouring, meaning it can help shape your body without the need for surgery. This makes it a safer and more comfortable option compared to traditional methods like liposuction.

## Triple Working Modes

- Simple mode: start treatment quickly and safely
- Expert mode: storable programs , helping professional beauticians and doctors memorize treatment plan
- AI mode: the clinical data collected by Omni over many years can be used directly by novices without worry
- 10 AI clinical protocols, 500 storeable customized treatment protocols
- True color touch screen, human-machine dialogue




## Specification




Model	RD_CET100
Frequency	448KHz
Treatment Handle	2pcs: CAP and RES
Treatment Tip	4pcs for CAP and 4pcs for RES
Display	8.4" colorful touch screen
Machine Weight	15kg




# Shanghai Omni Laser Skinology CO.,LTD

 Building 12-A-6F, Liando U Valley Songjiang International Enterprise Port Park, No. 1698 Minyi Road, Songjiang District, Shanghai, China

 Tel: +86-21-54847192 54847783

 [info@omni-laser.com](mailto:info@omni-laser.com)

 [www.omni-lasers.com](http://www.omni-lasers.com)

 [www.youtube.com/  
@OmniLaserSkinology](https://www.youtube.com/@OmniLaserSkinology)

